



TALENT
S U R G I C A L



CERVICAL SURGICAL INSTRUMENTS

We make more than 1385 surgical instruments.

We are proud to manufacture premium quality instruments.

Certification:



Registered with Therapeutic Goods Australia (Scissors)
Registered with Therapeutic Goods Australia (Needle Holder)
Registered with Therapeutic Goods Australia (Forceps)

www.tsurgicalinstruments.com



TALENT


S U R G I C A L




Instruments Name	Pages Number
DAVIS Vascular Spatula, Double Ended Dlight Curved b/b, 9 5/8 Inch	004
TOENNIS Vascular Spatula, Double Ended Tip Slight Curved, Tip size 3+4 m length 9 ½ Inch	
CUSHING Vessel Hook 1 x 4 BLUNT 7 ½ Inch	
KRAYENBUEHL Nerve Hook Tip SHORT BALL-POINT, 7 ¼ inch	005
LANDOLT-REULEN Dissector, Bayonet-Shaped 2.0 Diameter, Working Length 5 1/8 ,11 inch	
FAHLBUSCH Curette, Bayonet-Shaped 2.5 Diamter, MALLEABLE 6 1/8 working length, 10 ¼ Inch	
FERGUSSON Suction Tube, Luer Hub 2.0 Dia, Working 4 ¼, Total 7 ¾ Inch	006
FERGUSSON Suction Tube, Luer Hub 2.0 Dia, Working 5 1/8 , Total 8 ½ Inch	
FERGUSSON Suction Tube, Luer Hub 4.0 Dia, Working 6 1/4 , Total 9 5/8 Inch	
FERGUSSON Suction Tube, Luer Hub 3.0 Dia, Working 7 7/8 , Total 11 Inch	007
LOVE Nerve Root Retractor, Length 8 ¼ Inch	
AL-MEFTY Spoon, Tip STR, SHARP , Tip mm 2.8 Handle BLACK CAP, 7 7/8 Inch	
AL-MEFTY Spoon, Tip Angled, SHARP , Tip mm 2.8 Handle BLACK CAP, 7 7/8 Inch	008
KERRISON Bayonet-Shaped Bone Punch 40° upbiting. Jaw size 2.0, Jaw Opening 10 , 9 ½ Inch	
KERRISON Bayonet-Shaped Bone Punch, 40° upbiting. Jaw size 3.0, Jaw Opening 10 , 9 ½ Inch	
KERRISON Bayonet-Shaped Bone Punch 40° upbiting. Jaw size 4.0, Jaw Opening 10 , 9 ½ Inch	009
Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 1 mm, Working Length 7 Inch. Thin Model	
Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 1.5 mm, Working Length 7 Inch. Thin Model	
Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 2 mm, Working Length 7 Inch. Thin Model	010
Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 2.5 mm, Working Length 7 Inch. Thin Model	
Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 3 mm, Working Length 7 Inch. Thin Model	

Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 4 mm, Working Length 7 Inch. Thin Model	011
Black Coated Micro Pituitary Rongeur Jaw size 2 x 6, STR 7 inch working length	
Noir Detachable Rongeur, Straight Jaw, 1.5 mm 6 inch working length	
Noir Detachable Rongeur, Straight Jaw, 2 mm 9 inch working length	012
Noir Detachable Rongeur, Straight Jaw, 2 mm 9 inch working length	
RUSKIN Bone Rongeur, Double Action Jaw 2 mm, Straight mini, 6 ¼ Inch	
STELLBRINK Synovectomy Rongeur Jaw size 2.0 x 39.5 , Shape Jaw SLIGHT-CVD 6 7/8 inch	013
LEMPERT Bone Curette Tip 2.4 mm, 8 ½ Inch	
Bone Tamp, 2 Dia 6 ¼ inch	
Bone Tamp, 8 Dia 6 ¼ inch	014
Mallet, 155gr, 25 Dia, 7 ½ Inch	
LEXER-MINI Osteotome W/ PEEK Handle Tip 4 mm, size 7 Inch	
CLOWARD Style Cervical Retractor W/ Lip, width 18 mm, 7 inch	015
CLOWARD Style Cervical Retractor W/O Lip width 18 mm, 7 inch	
Kidney Tray, Length 6 ¾ Inch	
Metal Iodine Cup 175 ml	016
IODINE CUP, 400 ml	
Scalpel Handle #4, 8 7/8 Inch	
Scalpel Handle #4, 5 ¼ Inch	017
Scalpel Handle #3, 5 Inch	
WEITLANER Retractor, Blunt 3 x 4 Prong 5 1/8 Inch	
KOCHER-LANGENBECK Retractor Blade 70 x 14 mm, 8 ½ Inch	018
KOCHER-LANGENBECK Retractor Blade 35 x 15 mm, 8 ½ Inch	
BACKHAUS Towel Clamp 3 ½ Inch	
Non-Perforating Towel Forceps length 4 inch	019
KANTROWITZ Right Angle Forceps. 7 7/8 Inch	


JACOBSON Hemostatic Forceps Slight Curved 7 ½ Inch	
Mosquito Forceps Curved 6 Inch	020
COLLER Artery Forceps, Curved 6 ½ Inch	
LERICHE Forceps, Delicate Straight, 1 x 2 Teeth, 6 Inch	
KELLY Sponge Forceps, Jaw Curved, Without Ratchet, 12 ½ Inch	021
BABCOCK Forceps, 8 ½ inch	
Allis Atraumatic Forceps, Length 6 ¼ Inch	
COOLEY Atraumatic Vascular Tissue Forceps, Tip 2 mm, Length 7 ¾ Inch	022
ADSON Forceps, Delicate Smooth 2 x 3 teeth, Length 4 ¾ Inch	
JANSEN Bayonet Forceps 6 ¼ Inch	
Tissue Forceps, Standard Pattern, Teeth 1 x 2 Length 4 ½ inch	023
Dressing Forceps, Delicate, 5 1/8 inch	
TC MAYO-HEGAR Needle Holder, Straight 6 Inch	
DUROGRIP TC RYDER Needle Holder, smooth 7 inch	024
KAYE SuperCut Scissors, Delicate, Curved 7 ¾ inch. Black ring coated	
Noir METZENBAUM Dissecting Scissors B/B Curved 5 ¾ Inch	
MAYO-STILLE Dissecting Scissors, Beveled Blade, B/B Straight 5 ½ Inch	025
LISTER Bandage Scissors W/ One Large Ring, 8 Inch	
Steel Box	




- **DAVIS Vascular Spatula, Double Ended Dlight Curved b/b**
- Qty : 1
- Size : 9 5/8"
- Item Code : TCSI-01
- Japanese J1 ASTM F-899




- **TOENNIS Vascular Spatula, Double Ended Tip Slight Curved**
- Qty : 1
- Tip size 3+4 m length 9 1/2 Inch
- Item Code : TCSI-02
- Japanese J1 ASTM F-899




- **CUSHING Vessel Hook 1 x 4 BLUNT 7 1/2 Inch**
- Qty : 1
- Item Code : TCSI-03
- Japanese J1 ASTM F-899




- **KRAYENBUEHL Nerve Hook**
Tip SHORT BALL-POINT
- Qty : 1
- Size : 7 1/4"
- Item Code : TCSI-04
- Japanese J1 ASTM F-899




- **LANDOLT-REULEN Dissector, Bayonet-Shaped**
2.0 Diameter, Working Length 5 1/8 ,11 inch
- Qty : 1
- Size : 7"
- Item Code : TCSI-05
- Japanese J1 ASTM F-899




- **FAHLBUSCH Curette, Bayonet-Shaped**
2.5 Diameter, MALLEABLE 6 1/8 working length, 10 1/4 Inch
- Qty : 1
- Item Code : TCSI-06
- Japanese J1 ASTM F-899




- **FERGUSSON Suction Tube, Luer Hub
2.0 Dia, Working 4 1/4, Total 7 3/4 Inch**
- Qty : 1
- Item Code : TCSI-07
- Japanese J1 ASTM F-899




- **FERGUSSON Suction Tube, Luer Hub
2.0 Dia, Working 5 1/8, Total 8 1/2 Inch**
- Qty : 1
- Item Code : TCSI-08
- Japanese J1 ASTM F-899




- **FERGUSSON Suction Tube, Luer Hub
4.0 Dia, Working 6 1/4, Total 9 5/8 Inch**
- Qty : 1
- Item Code : TCSI-09
- Japanese J1 ASTM F-899



- **FERGUSSON Suction Tube, Luer Hub
3.0 Dia, Working 7 7/8 , Total 11 Inch**
- Qty : 1
- Item Code : TCSI-10
- Japanese J1 ASTM F-899




- **LOVE Nerve Root Retractor,
Length 8 1/4 Inch**
- Qty : 1
- Item Code : TCSI-11
- Japanese J1 ASTM F-899




- **AL-MEFTY Spoon,
Tip STR, SHARP , Tip mm 2.8 Handle
BLACK CAP, 7 7/8 Inch**
- Qty : 1
- Item Code : TCSI-12
- Japanese J1 ASTM F-899




- **AL-MEFTY Spoon, Tip Angled, SHARP , Tip mm 2.8 Handle BLACK CAP, 7 7/8 Inch**
- Qty : 1
- Item Code : TCSI-13
- Japanese J1 ASTM F-899




- **KERRISON Bayonet-Shaped Bone Punch 40° upbiting. Jaw size 2.0, Jaw Opening 10 , 9 ½ Inch**
- Qty : 1
- Item Code : TCSI-14
- Japanese J1 ASTM F-899




- **KERRISON Bayonet-Shaped Bone Punch 40° upbiting. Jaw size 3.0, Jaw Opening 10 , 9 ½ Inch**
- Qty : 1
- Item Code : TCSI-15
- Japanese J1 ASTM F-899



- **KERRISON Bayonet-Shaped Bone Punch**
40° upbiting. Jaw size 4.0, Jaw Opening 10 , 9 ½ Inch
- Qty : 1
- Item Code : TCSI-16
- Japanese J1 ASTM F-899



- **Kerrison Bone Punches, Non Detachable, Silver.**
Jaw Size 1 mm, Working Length 7 Inch. Thin Model
- Qty : 1
- Item Code : TCSI-17
- Japanese J1 ASTM F-899



- **Kerrison Bone Punches, Non Detachable, Silver.**
Jaw Size 1.5 mm, Working Length 7 Inch. Thin Model
- Qty : 1
- Item Code : TCSI-18
- Japanese J1 ASTM F-899



- **Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 2 mm, Working Length 7 Inch. Thin Model**
- Qty : 1
- Item Code : TCSI-19
- Japanese JI ASTM F-899




- **Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 2.5 mm, Working Length 7 Inch. Thin Model**
- Qty : 1
- Item Code : TCSI-20
- Japanese JI ASTM F-899



- **Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 3 mm, Working Length 7 Inch. Thin Model**
- Qty : 1
- Item Code : TCSI-21
- Japanese JI ASTM F-899



- **Kerrison Bone Punches, Non Detachable, Silver. Jaw Size 4 mm, Working Length 7 Inch. Thin Model**
- Qty : 1
- Item Code : TCSI-22
- Japanese J1 ASTM F-899



- **Black Coated Micro Pituitary Rongeur**
- **Jaw size 2 x 6, STR 7 inch working length**
- Qty : 1
- Item Code : TCSI-23
- Japanese J1 ASTM F-899



- **Noir Detachable Rongeur**
- **Straight Jaw, 1.5 mm**
- **6 inch working length**
- Qty : 1
- Item Code : TCSI-24
- Japanese J1 ASTM F-899



- **Noir Detachable Rongeur
Straight Jaw, 2 mm
9 inch working length**
- Qty : 1
- Item Code : TCSI-25
- Japanese J1 ASTM F-899




- **Noir Detachable Rongeur
Straight Jaw, 2 mm
9 inch working length**
- Qty : 1
- Item Code : TCSI-26
- Japanese J1 ASTM F-899




- **RUSKIN Bone Rongeur, Double Action
Jaw 2 mm, Straight mini, 6 1/4 Inch**
- Qty : 1
- Item Code : TCSI-27
- Japanese J1 ASTM F-899




- **STELLBRINK Synovectomy Rongeur**
Jaw size 2.0 x 39.5 , Shape Jaw SLIGHT-CVD 6 7/8 inch
- Qty : 1
- Item Code : TCSI-28
- Japanese J1 ASTM F-899



- **LEMPERT Bone Curette**
Tip 2.4 mm, 8 ½ Inch
- Qty : 1
- Item Code : TCSI-29
- Japanese J1 ASTM F-899




- **Bone Tamp**
2 Dia 6 ¼ inch
- Qty : 1
- Item Code : TCSI-30
- Japanese J1 ASTM F-899



- **Bone Tamp**
8 Dia 6 ¼ inch
- Qty : 1
- Item Code : TCSI-31
- Japanese J1 ASTM F-899



- **Mallet**
155gr, 25 Dia, 7 ½ Inch
- Qty : 1
- Item Code : TCSI-32
- Japanese J1 ASTM F-899



- **LEXER-MINI Osteotome W/ PEEK Handle**
Tip 4 mm, size 7 Inch
- Qty : 1
- Item Code : TCSI-33
- Japanese J1 ASTM F-899



- **CLOWARD Style Cervical Retractor W/ Lip , width 18 mm, 7 inch**
- Qty : 1
- Item Code : TCSI-34
- Japanese J1 ASTM F-899



- **CLOWARD Style Cervical Retractor W/O Lip width 18 mm, 7 inch**
- Qty : 1
- Item Code : TCSI-35
- Japanese J1 ASTM F-899



- **Kidney Tray**
- **Length 6 ¾ Inch**
- Qty : 1
- Item Code : TCSI-36
- Japanese J1 ASTM F-899



- **Metal Iodine Cup 175 ml**
- Qty : 1
- Item Code : TCSI-37
- Japanese J1 ASTM F-899

S U R G I C A L



- **IODINE CUP, 400 ml**
- Qty : 1
- Item Code : TCSI-38
- Japanese J1 ASTM F-899

S U R G I C A L




- **Scalpel Handle #4
8 7/8 Inch**
- Qty : 1
- Item Code : TCSI-39
- Japanese J1 ASTM F-899


S U R G I C A L




- **Scalpel Handle #4**
5 1/4 Inch
- Qty : 1
- Item Code : TCSI-40
- Japanese J1 ASTM F-899




- **Scalpel Handle #3**
5 Inch
- Qty : 1
- Item Code : TCSI-41
- Japanese J1 ASTM F-899




- **WEITLANER Retractor, Blunt 3 x 4 Prong**
5 1/8 Inch
- Qty : 2
- Item Code : TCSI-42
- Japanese J1 ASTM F-899




- **KOCHER-LANGENBECK Retractor**
Blade 70 x 14 mm, 8 ½ Inch
- Qty : 2
- Item Code : TCSI-43
- Japanese J1 ASTM F-899




- **KOCHER-LANGENBECK Retractor**
Blade 35 x 15 mm, 8 ½ Inch
- Qty : 2
- Item Code : TCSI-44
- Japanese J1 ASTM F-899




- **BACKHAUS Towel Clamp**
3 ½ Inch
- Qty : 4
- Item Code : TCSI-45
- Japanese J1 ASTM F-899




- **Non-Perforating Towel Forceps**
length 4 inch
- Qty : 1
- Item Code : TCSI-46
- Japanese J1 ASTM F-899




- **KANTROWITZ Right Angle Forceps**
7 7/8 Inch
- Qty : 1
- Item Code : TCSI-47
- Japanese J1 ASTM F-899



- **JACOBSON Hemostatic Forceps**
Slight Curved 7 1/2 Inch
- Qty : 2
- Item Code : TCSI-48
- Japanese J1 ASTM F-899




- **Mosquito Forceps Curved 6 Inch**
- Qty : 4
- Item Code : TCSI-49
- Japanese J1 ASTM F-899




- **COLLER Artery Forceps, Curved 6 ½ Inch**
- Qty : 4
- Item Code : TCSI-50
- Japanese J1 ASTM F-899




- **LERICHE Forceps, Delicate Straight 1 x 2 Teeth, 6 Inch**
- Qty : 2
- Item Code : TCSI-51
- Japanese J1 ASTM F-899



- **KELLY Sponge Forceps, Jaw Curved, Without Ratchet, 12 ½ Inch**
- Qty : 2
- Item Code : TCSI-52
- Japanese J1 ASTM F-899




- **BABCOCK Forceps**
- **8 ½ inch**
- Qty : 2
- Item Code : TCSI-53
- Japanese J1 ASTM F-899




- **Allis Atraumatic Forceps**
- **Length 6 ¼ Inch**
- Qty : 2
- Item Code : TCSI-54
- Japanese J1 ASTM F-899




- **COOLEY Atraumatic Vascular Tissue Forceps**
Tip 2 mm, Length 7 ¾ Inch
- Qty : 1
- Item Code : TCSI-55
- Japanese J1 ASTM F-899




- **ADSON Forceps**
Delicate Smooth 2 x 3 teeth, Length 4 ¾ Inch
- Qty : 1
- Item Code : TCSI-56
- Japanese J1 ASTM F-899



- **JANSEN Bayonet Forceps**
6 ¼ Inch
- Qty : 1
- Item Code : TCSI-57
- Japanese J1 ASTM F-899



- **Tissue Forceps, Standard Pattern
Teeth 1 x 2 Length 4 1/2 inch**
- Qty : 1
- Item Code : TCSI-58
- Japanese J1 ASTM F-899




- **Dressing Forceps, Delicate
5 1/8 inch**
- Qty : 2
- Item Code : TCSI-59
- Japanese J1 ASTM F-899



- **TC MAYO-HEGAR Needle Holder, Straight
6 Inch**
- Qty : 1
- Item Code : TCSI-60
- Japanese J1 ASTM F-899



- **CDUROGRIP TC RYDER Needle Holder, smooth 7 inch**
- Qty : 1
- Item Code : TCSI-61
- Japanese J1 ASTM F-899




- **KAYE SuperCut Scissors, Delicate, Curved 7 3/4 inch. Black ring coated**
- Qty : 1
- Item Code : TCSI-62
- Japanese J1 ASTM F-899



- **Noir METZENBAUM Dissecting Scissors B/B Curved 5 3/4 Inch**
- Qty : 1
- Item Code : TCSI-63
- Japanese J1 ASTM F-899



- **MAYO-STILLE Dissecting Scissors**
Beveled Blade, B/B Straight 5 ½ Inch
- Qty : 1
- Item Code : TCSI-64
- Japanese J1 ASTM F-899



- **LISTER Bandage Scissors**
W/ One Large Ring, 8 Inch
- Qty : 1
- Item Code : TCSI-65
- Japanese J1 ASTM F-899



- **Steel Box**
- Qty : 1
- Special box TCSI-66
- Japanese J1 ASTM F-899



Sialkot, Pakistan Office

Contact Person: Hasham khalid



Address: A/1, Small Industries Estate, Near SIE Telephone Exchange PTCL, Sialkot – Punjab 51310, Pakistan.



Yakshan Khalid: +92 331 478 6083



info@tsurgicalinstruments.com
Talentsurgicalaustralia@gmail.com



Australia Office:

Contact Person: Hasham khalid



Pelham Crescent, Wyndham Vale, VIC 3024 Australia



Hasham Khalid: +61 430 855 424 / +61 481 335 786



info@tsurgicalinstruments.com
Talentsurgicalaustralia@gmail.com



TALENT
S U R G I C A L

